

# Mi-OPA<sup>®</sup>

## High Level Disinfectant for Medical Devices




- Rapid action.
- Ready to use.
- Non rusting.
- Excellent material compatibility.
- Flexible reprocessing.
- Proven Safety.
- Long Life.

### Active Ingredients

Ortho-phthalaldehyde 0.55% w/v  
with chelating agents, buffer, corrosion  
inhibitors & aqua



**FMS**  
You Say. We Disinfect.



## High Level Disinfectant for Medical Devices

### AREAS OF APPLICATION:

Mi-OPA is a high level disinfectant intended to disinfect heat sensitive reusable medical devices. It has excellent material compatibility and is suitable for reprocessing of instruments made of rubber, plastic metal, glass and porcelain. It can be used for manual, semi-automatic as well as completely automatic procedures. It can also be used in ultrasonic baths.

### DIRECTIONS FOR USE :

**Manual Processing :** immerse device completely in Mi-OPA filling all lumens and eliminating air pockets Immerse for minimum 5 minutes to achieve high level disinfection.

**Re-processing using Automatic Endoscope Re-processor:** Set the temperature to 25°C with an immersion time of 5 minutes.



After removing the device from Mi-OPA, wash it thoroughly with large volume of sterile water.

### PRECAUTIONS :

- Store at controlled room temperature, 15° to 30°C.
- For external use only.
- Keep soaking trays covered at all times.
- Avoid direct contact with skin, eyes and clothing.
- Wear personal protective gear like rubber gloves, eye shield while handling.
- Wash hands thoroughly with soap and water after handling.
- Flush skin, eyes with fresh water in case of contact.
- In case of emergency seek medical help.
- Keep out of reach of children.

**The solution in secondary container can be used for a period of 14 days from the date of transfer.**

Pack Size : 5 Liter



You Say. We Disinfect.

### FAITH MICROSOLUTIONS

ISO 9001 : 2015 & ISO 13485 : 2016 Company

**Head office:** 27, Ashoka Chambers, 5-B Rajendra Park,  
New Delhi-110060, INDIA

**Phone :** 011-47479797 (20 Lines)

**Customer Care :** info@faithmicrosolutions.com

**Url :** www.faithmicrosolutions.com

Authorised Representative :



Connect with us at :